

FOXRED1 antibody



Catalog Number: 110716

Product name

FOXRED1 antibody

Immunogen

[Human FOXRED1 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

FOXRED1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 30 kDa, 50 kDa

Preparation

This antibody was obtained by immunization of FOXRED1 recombinant protein (Accession Number: NM_017547). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110716(FOXRED1 Antibody) at dilution of 1:300