

MNK1 antibody



Catalog Number: 112716

Product name

MNK1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MNK1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, HepG2 cells, HL-60 cells, mouse small intestine tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 51 kDa, 47 kDa

Preparation

This antibody was obtained by immunization of MNK1 recombinant protein (Accession Number: BC002755). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112716(MKNK1 antibody) at dilution of 1:600



Immunofluorescent analysis of HepG2 cells using Catalog No:112716(MKNK1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)