

Phospho-MAP3K5-S83 pAb



Catalog Number: 166426

Product name

Phospho-MAP3K5-S83 pAb

Specificity

Human

Antibody description

Polyclonal antibody to Phospho-MAP3K5-S83

Preparation

Antigen: A phospho specific peptide corresponding to residues surrounding S83 of human MAP3K5

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Validations



Western blot - Phospho-MAP3K5-S83 pAb

Western blot analysis of extracts from K562 cells using Phospho-MAP3K5-S83 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% BSA.



Immunohistochemistry - Phospho-MAP3K5-S83 pAb

Immunohistochemistry of paraffin-embedded human breast carcinoma using Phospho-MAP3K5-S83 antibody .