

NUDT15 Polyclonal Antibody



Catalog Number: 165826

Product name

NUDT15 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to NUDT15

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-164 of human NUDT15 (NP_060753.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Validations




Western blot - NUDT15 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using NUDT15 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - NUDT15 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver using NUDT15 antibody at dilution of



1:100 (40x lens).

NUDT15 Polyclonal Antibody



Catalog Number: 165826

Immunohistochemistry - NUDT15 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using NUDT15 antibody at dilution of 1:100 (40x lens).