

SERCA2 antibody



Catalog Number: 115126

Product name

SERCA2 antibody

Immunogen

[Human SERCA2,ATP2A2 Recombinant protein \(His tag\)](#)

Specificity

Human,Mouse,Rat; other species not tested.

Antibody description

SERCA2 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 100kd

Preparation

This antibody was obtained by immunization of SERCA2 recombinant protein (Accession Number: NM_170665). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:1000-1:10000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115126(SERCA2,ATP2A2 antibody) at dilution of 1:3000