

Anti-BRK1 Rabbit antibody



Catalog Number: 167526

Product name

Anti-BRK1 Rabbit antibody

Specificity

Human

Antibody description

BRK1 Rabbit polyclonal antibody

Preparation

Antigen: Synthetic peptide of human BRK1

Formulation

-20°C, pH7.4 PBS, 0.05% NaN₃, 40% Glycerol

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC-P

Dilutions


WB: 1/500-1/2000

IHC: 1/50-1/200

Validations



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 167526(BRK1 Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)



The image is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using 167526(BRK1 Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original



magnification: ×200)

Gel: 10%SDS-PAGE, Lysate: 40 µg, Lane: A172 cells, Primary antibody: 167526(BRK1 Antibody) at dilution 1/350, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1

Anti-BRK1 Rabbit antibody



Catalog Number: 167526

minute