

# RRN3 antibody



Catalog Number: 114926

## Product name

RRN3 antibody

## Immunogen

[Human RRN3 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RRN3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells, K-562 cells, mouse testis tissue. Observed molecular weight by Western-blot: 74 kDa

## Preparation

This antibody was obtained by immunization of RRN3 recombinant protein (Accession Number: NM\_018427). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114926(RRN3 Antibody) at dilution of 1:600