

# EHD4 antibody



Catalog Number: 110226

## Product name

EHD4 antibody

## Immunogen

[Human EHD4 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

EHD4 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse testis tissue, HeLa cells. Observed molecular weight by Western-blot: 61 kDa

## Preparation

This antibody was obtained by immunization of EHD4 recombinant protein (Accession Number: NM\_139265). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110226(EHD4 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using EHD4 antibody Catalog No:110226 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).