

MIP antibody



Catalog Number: 112626

Product name

MIP antibody

Immunogen

[Human MIP Recombinant protein \(GST tag & His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

MIP Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 26-28kd

Preparation

This antibody was obtained by immunization of MIP recombinant protein (Accession Number: NM_012064). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:112626(MIP antibody) at dilution of 1:500