

ATRAID antibody



Catalog Number: 108326

Product name

Polyclonal

ATRAID antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

ATRAID Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HEK-293 cells, mouse kidney tissue. Observed molecular weight by Western-blot: 28-35 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of ATRAID recombinant protein (Accession Number: NM_016085). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108326(C2orf28 Antibody) at dilution of 1:300