

THADA antibody



Catalog Number: 116036

Product name

THADA antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

THADA Rabbit Polyclonal antibody. Positive WB detected in A549 cells, CHO cells, HeLa cells. Positive IHC detected in human pancreas tissue, human testis tissue, mouse ovary tissue, rat ovary tissue. Observed molecular weight by Western-blot: 183kd

Preparation

This antibody was obtained by immunization of THADA recombinant protein (Accession Number: BC025773). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations




WB result of Catalog No:116036 from Takashi



Kobayashi.

A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:116036(THADA antibody) at dilution of 1:300




Immunohistochemical of paraffin-embedded human pancreas using Catalog No:116036(THADA antibody) at dilution of 1:50 (under 10x lens)

THADA antibody



Catalog Number: 116036



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:116036(THADA antibody) at dilution of 1:50 (under 40x lens)