

# RAB38 antibody



Catalog Number: 114436

## Product name

RAB38 antibody

## Immunogen

[Human RAB38 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

human, mouse; other species not tested.

## Antibody description

RAB38 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 27 kDa

## Preparation

This antibody was obtained by immunization of RAB38 recombinant protein (Accession Number: NM\_022337). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114436(RAB38 antibody) at dilution of 1:500