

Anti-APPL Rabbit antibody



Catalog Number: 167436

Product name

Anti-APPL Rabbit antibody

WB: 1/1000-1/5000

IHC: 1/20

Specificity

Human, Mouse, Rat

ICC/IF: 1/50

Antibody description

APPL Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human APPL

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF

Dilutions

Validations



Western blot detection of APPL in Rat Brain, C6, 3T3, Hela cell lysates using APPL Rabbit pAb (1:1000 diluted). Predicted band size: 80 kDa. Observed band size: 80 kDa.



Immunohistochemistry of APPL in paraffin-embedded Human lung cancer tissue using APPL Rabbit pAb at dilution 1/20