

SF3B14 antibody



Catalog Number: 115136

Product name

Polyclonal

SF3B14 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

SF3B14 Rabbit Polyclonal antibody. Positive WB detected in mouse uterus tissue. Observed molecular weight by Western-blot: 14kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of SF3B14 recombinant protein (Accession Number: NM_016047). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

mouse uterus tissue were subjected to SDS PAGE followed by western blot with Catalog No:115136(SF3B14 antibody) at dilution of 1:300

