

PDRG1 antibody



Catalog Number: 113636

Product name

PDRG1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PDRG1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Positive IF detected in HepG2 cells. Positive IHC detected in human kidney tissue, human liver tissue. Observed molecular weight by Western-blot: 16 kDa

Preparation

This antibody was obtained by immunization of PDRG1 recombinant protein (Accession Number: NM_030815). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113636(PDRG1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human kidney using Catalog No:113636(PDRG1 antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded human kidney using Catalog No:113636(PDRG1

PDRG1 antibody



Catalog Number: 113636

antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using PDRG1 antibody Catalog No:113636 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).