

# THOC1 antibody



Catalog Number: 116046

## Product name

THOC1 antibody

## Immunogen

[Human THOC1 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

THOC1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 84 kDa

## Preparation

This antibody was obtained by immunization of THOC1 recombinant protein (Accession Number: NM\_005131). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116046(THoc1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116046(THoc1 antibody) at dilution of 1:100 (under 10x lens)