

TALPID3 antibody



Catalog Number: 115846

Product name

TALPID3 antibody

Immunogen

[Human KIAA0586 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

TALPID3 Rabbit Polyclonal antibody. Positive WB detected in mouse embryo tissue. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of TALPID3 recombinant protein (Accession Number: NM_014749). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



mouse embryo tissue were subjected to SDS PAGE followed by western blot with Catalog No:115846(KIAA0586 Antibody) at dilution of 1:600