

SCNN1A Polyclonal Antibody



Catalog Number: 164546

Product name

SCNN1A Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to SCNN1A

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 105-225 of human SCNN1A (NP_001029.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:200 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - SCNN1A Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SCNN1A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.





Immunofluorescence - SCNN1A Polyclonal Antibody

Immunofluorescence analysis of human lung cancer cells using SCNN1A antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

SCNN1A Polyclonal Antibody




Catalog Number: 164546



Immunofluorescence - SCNN1A Polyclonal Antibody

Immunofluorescence analysis of rat lung cells using SCNN1A antibody at dilution of 1:100. Blue: DAPI for nuclear staining.




DAPI for nuclear staining.

Immunofluorescence - SCNN1A Polyclonal Antibody

Immunofluorescence analysis of mouse kidney cells using SCNN1A antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Immunofluorescence - SCNN1A Polyclonal Antibody

Immunofluorescence analysis of human lung cancer cells using SCNN1A antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence - SCNN1A Polyclonal Antibody

Immunofluorescence analysis of rat lung cells using SCNN1A antibody at dilution of 1:100. Blue: DAPI for nuclear staining.