

NAP1L4 antibody



Catalog Number: 112946

Product name

NAP1L4 antibody

Recommended Dilution:

WB: 1:1000-1:10000

Immunogen

IP: 1:500-1:5000

[Human NAP1L4 Recombinant protein \(His tag\)](#)

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

NAP1L4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, mouse liver tissue, mouse testis tissue. Positive IP detected in HeLa cells. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 52 kDa



Preparation

This antibody was obtained by immunization of NAP1L4 recombinant protein (Accession Number: NM_001369376). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112946(NAP1L4 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP


Dilutions




Immunohistochemical of paraffin-embedded human liver using Catalog No:112946(NAP1L4 antibody) at dilution of 1:100 (under 10x lens)

NAP1L4 antibody

Catalog Number: 112946



Immunohistochemical of paraffin-embedded human liver using Catalog No:112946(NAP1L4 antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-NAP1L4 (IP:Catalog No:112946, 4ug; Detection:Catalog No:112946 1:1000) with HeLa cells lysate 4000ug.