

SIRT5 antibody



Catalog Number: 115246

Product name

SIRT5 antibody

Immunogen

[Human SIRT5 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SIRT5 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse liver tissue. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 33kd

Preparation

This antibody was obtained by immunization of SIRT5 recombinant protein (Accession Number: NM_001376801). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115246(SIRT5 antibody) at dilution of 1:400



IP Result of anti-SIRT5 (IP:Catalog No:115246, 4ug; Detection:Catalog No:115246 1:500) with mouse liver tissue lysate 6000ug.