

CEBPB antibody



Catalog Number: 109246

Product name

CEBPB antibody

Immunogen

[Human CEBPB Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Rat; other species not tested.

Antibody description

CEBPB Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells. Observed molecular weight by Western-blot: 36 kDa

Preparation

This antibody was obtained by immunization of CEBPB recombinant protein (Accession Number: NM_005194). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:109246(CEBPB Antibody) at dilution of 1:300