

ATP1B3 antibody



Catalog Number: 108346

Product name

ATP1B3 antibody

Specificity

human; other species not tested.

Antibody description

ATP1B3 Rabbit Polyclonal antibody. Positive WB detected in human lung tissue. Positive IHC detected in human prostate cancer. Observed molecular weight by Western-blot: 60-64 kDa

Preparation

This antibody was obtained by immunization of ATP1B3 recombinant protein (Accession Number: NM_001679). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:108346(ATP1B3 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:108346(ATP1B3 antibody) at dilution of 1:50 (under 10x lens)