

ACOT6 antibody



Catalog Number: 107746

Product name

ACOT6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ACOT6 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HeLa cells. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 23kd

Preparation

This antibody was obtained by immunization of ACOT6 recombinant protein (Accession Number: NM_001365789). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107746(ACOT6 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human kidney slide using Catalog No:107746(ACOT6 Antibody) at dilution of 1:50