

Anti-PDPK1 Rabbit antibody



Catalog Number: 167056

Product name

Anti-PDPK1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

PDPK1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human PDPK1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions


WB: 1/1000

IHC: 1/20

Validations



Immunohistochemistry of PDPK1 in paraffin-embedded Human colon cancer tissue using



PDPK1 Rabbit pAb at dilution 1/20

Western blot detection of PDPK1 in HeLa cell lysates using PDPK1 Rabbit pAb(1:1000 diluted). Predicted band size:63kDa.Observed band size:63kDa.