

SIN3A Polyclonal Antibody



Catalog Number: 161856

Product name

SIN3A Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to SIN3A

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human SIN3A (NP_056292.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - SIN3A Polyclonal Antibody

Western blot analysis of extracts of 293T cells, using SIN3A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - SIN3A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat liver using SIN3A antibody at dilution of 1:100



(40x lens).


SIN3A Polyclonal Antibody



Catalog Number: 161856

Immunohistochemistry - SIN3A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using SIN3A antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - SIN3A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using SIN3A antibody at dilution of 1:100 (40x lens).