

DAK antibody



Catalog Number: 109856

Product name

DAK antibody

Immunogen

[Human DAK Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DAK Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, mouse uterus tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 59kd

Preparation

This antibody was obtained by immunization of DAK recombinant protein (Accession Number: XM_017017521). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109856(DAK antibody) at dilution of 1:500



IP Result of anti-DAK (IP:Catalog No:109856, 3ug; Detection:Catalog No:109856 1:500) with HeLa cells lysate 2150ug.