

TRIP12 Polyclonal Antibody



Catalog Number: 166256

Product name

TRIP12 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to TRIP12

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human TRIP12 (NP_001271143.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:1000

IHC 1:50 - 1:100


IF 1:50 - 1:100

Validations



Western blot - TRIP12 Polyclonal Antibody


Western blot analysis of extracts of various cell lines, using TRIP12 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 30s.

Immunohistochemistry - TRIP12 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon using TRIP12 antibody at dilution of



1:100 (40x lens).


TRIP12 Polyclonal Antibody



Catalog Number: 166256

Immunohistochemistry - TRIP12 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat testis using TRIP12 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - TRIP12 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse testis using TRIP12 antibody at dilution of 1:100 (40x lens).