

ARG2 antibody



Catalog Number: 108256

Product name

ARG2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ARG2 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue, human kidney tissue. Positive WB detected in HEK-293 cells, mouse liver tissue, rat liver tissue. Observed molecular weight by Western-blot: 39-42kd

Preparation

This antibody was obtained by immunization of ARG2 recombinant protein (Accession Number: NM_001172). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108256(ARG2 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:108256(ARG2 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:108256(ARG2 antibody) at dilution of 1:100

ARG2 antibody



Catalog Number: 108256

(under 40x lens)