

Fibrinogen beta chain antibody



Catalog Number: 110656

Product name

Fibrinogen beta chain antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Fibrinogen beta chain Rabbit Polyclonal antibody.
Positive IHC detected in human kidney tissue.
Positive WB detected in mouse heart tissue,
mouse brain tissue, mouse kidney tissue, mouse
liver tissue. Observed molecular weight by
Western-blot: 47-56 kDa

Preparation

This antibody was obtained by immunization of
Fibrinogen beta chain recombinant protein
(Accession Number: NM_005141). Purification
method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH
7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



mouse heart tissue were subjected to SDS PAGE
followed by western blot with Catalog
No:110656(FGB antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded
human kidney using Catalog No:110656(FGB
antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded
human kidney using Catalog No:110656(FGB
antibody) at dilution of 1:50 (under 40x lens)