

FXYD4 antibody



Catalog Number: 110756

Product name

FXYD4 antibody

IHC: 1:20-1:200


Specificity

Human, Mouse; other species not tested.

Validations

Antibody description

FXYD4 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, COLO 320 cells, human kidney tissue, human testis tissue, mouse kidney tissue. Positive IHC detected in mouse testis tissue. Observed molecular weight by Western-blot: 9-10 kDa




Preparation

This antibody was obtained by immunization of FXYD4 recombinant protein (Accession Number: NM_001184963). Purification method: Antigen affinity purified.

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110756(FXYD4 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunohistochemistry of paraffin-embedded mouse testis tissue slide using Catalog No:110756(FXYD4 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded mouse testis tissue slide using Catalog

FXYD4 antibody



Catalog Number: 110756

No:110756(FXYD4 Antibody) at dilution of 1:50 (under 40x lens)