

ZMAT2 antibody



Catalog Number: 117066

Product name

ZMAT2 antibody

Immunogen

[Human ZMAT2 Recombinant protein \(His tag\)](#)

Specificity

Human, Rat; other species not tested.

Antibody description

ZMAT2 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, C6 cells, HeLa cells, Jurkat cells, MCF-7 cells. Observed molecular weight by Western-blot: 24 kDa

Preparation

This antibody was obtained by immunization of ZMAT2 recombinant protein (Accession Number: NM_144723). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:117066(ZMAT2 Antibody) at dilution of 1:600