

ANKRD13A antibody



Catalog Number: 108066

Product name

ANKRD13A antibody

Immunogen

[Human ANKRD13A Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ANKRD13A Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 67 kDa

Preparation

This antibody was obtained by immunization of ANKRD13A recombinant protein (Accession Number: NM_033121). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108066(ANKRD13A Antibody) at dilution of



1:300

Immunofluorescent analysis of HepG2 cells using Catalog No:108066(ANKRD13A Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)