

TRPV3 Polyclonal Antibody



Catalog Number: 161166

Product name

TRPV3 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to TRPV3

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 400-500 of human TRPV3 (NP_659505.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Immunohistochemistry - TRPV3 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat heart using TRPV3 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - TRPV3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using TRPV3 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - TRPV3 Polyclonal Antibody

TRPV3 Polyclonal Antibody



Catalog Number: 161166

Immunohistochemistry of paraffin-embedded

mouse brain using TRPV3 antibody at dilution of 1:100 (40x lens).