

ZNF253 antibody



Catalog Number: 116966

Product name

ZNF253 antibody

Immunogen

[Human ZNF253 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ZNF253 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, SMMC-7721 cells. Observed molecular weight by Western-blot: 58-65 kDa

Preparation

This antibody was obtained by immunization of ZNF253 recombinant protein (Accession Number: NM_021047). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116966(ZNF253 Antibody) at dilution of 1:300