

DNMT1 antibody



Catalog Number: 109966

Product name

DNMT1 antibody

Immunogen

[Human DNMT1 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

DNMT1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 180-200 kDa

Preparation

This antibody was obtained by immunization of DNMT1 recombinant protein (Accession Number: NM_001130823). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:1000


IP: 1:200-1:2000

IF: 1:10-1:100

Validations




IP Result of anti-DNMT1 (IP:Catalog No:109966, 5ug; Detection:Catalog No:109966 1:500) with



HEK-293 cells lysate 2500ug.

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109966(DNMT1 Antibody) at dilution of 1:300



DNMT1 antibody



Catalog Number: 109966

Immunofluorescent analysis of HepG2 cells using

Catalog No:109966(DNMT1 Antibody) at dilution
of 1:25 and Rhodamine-Goat anti-Rabbit IgG