

# IFT25 antibody



Catalog Number: 111666

## Product name

IFT25 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

IFT25 Rabbit Polyclonal antibody. Positive IF detected in hTERT-RPE1 cells. Positive WB detected in HeLa cells, A549 cells, mouse brain tissue. Observed molecular weight by Western-blot: 25 kDa

## Preparation

This antibody was obtained by immunization of IFT25 recombinant protein (Accession Number: NM\_016126). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111666(HSPB11 antibody) at dilution of 1:500



IF result (the base of cilia stain) of anti-HSPB11 (IFT25, Catalog No:111666, 1:50) with serum-starved hTERT-RPE1 (PFA fixed) by Dr. Moshe Kim.