

# Anti-EGF (7C11) Mouse antibody



Catalog Number: 168076

## Product name

Monoclonal

Anti-EGF (7C11) Mouse antibody

## Ig Type

## Specificity

Mouse IgG1

Human

## Applications

## Antibody description

ELISA

EGF (7C11) Mouse monoclonal antibody

## Validations

## Preparation

Antigen: Purified recombinant Human EGF protein fragments expressed in E.coli.


## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality



Standard Curve for Human EGF: Capture Antibody Mouse mAb 201165-5E9 (5E9-A10-A6) to Human EGF at 4u03bcg/ml and Detector Antibody Mouse mAb 168076(7C11-G10-H8) to Human EGF at 1u03bcg/ml.