

Anti-SGSM1/RUTBC2 Rabbit antibody



Catalog Number: 167176

Product name

Anti-SGSM1/RUTBC2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

SGSM1/RUTBC2 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human SGSM1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/1000


ICC/IF: 1/20-1/100

FC: 1/20

Validations



Immunocytochemistry of SGSM1(green) in HeLa cells using SGSM1 Rabbit pAb at dilution 1/50, and DAPI(blue)



Western blot detection of SGSM1 in K562, HeLa cell lysates using SGSM1 Rabbit pAb(1:1000 diluted). Predicted band size: 130kDa. Observed band size: 130kDa.