

CREST antibody



Catalog Number: 107176

Product name

CREST antibody

Immunogen

[Human CREST Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

CREST Mouse Monoclonal antibody. Positive WB detected in MCF-7 cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of CREST recombinant protein (Accession Number: NM_015558). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107176(CREST Antibody) at dilution of 1:500