

Anti-BAP31 Rabbit antibody



Catalog Number: 168576

Product name

Anti-BAP31 Rabbit antibody

Dilutions

1:1000-1:5000

Specificity

1:50-1:200

Human, Mouse

1:20-1:50

Antibody description

1:20-1:100

BAP31 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human BAP31

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF, FC

Validations



Immunocytochemistry of BAP31(green) in HeLa cells using BAP31 Rabbit mAb at dilution 1/200,



and DAPI(blue)

Western blot detection of BAP31 in K562, HeLa cell lysates using BAP31 Rabbit mAb(1:1000 diluted). Predicted band size: 28kDa. Observed band size: 28kDa.