

BRCA1 antibody

Catalog Number: 108576

Product name

BRCA1 antibody

Immunogen

[Human BRCA1 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

BRCA1 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human lung cancer tissue. Positive WB detected in HeLa cells, A431 cells, HEK-293 cells, MCF-7 cells. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of BRCA1 recombinant protein (Accession Number: BC115037). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:117223(BRCA1 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:117223(BRCA1 Antibody) at dilution of 1:50