

GINS3 antibody



Catalog Number: 110976

Product name

Polyclonal

GINS3 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

GINS3 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 24kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of GINS3 recombinant protein (Accession Number: NM_022770). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110976(GINS3 antibody) at dilution of 1:800