

# CYFIP2 antibody



Catalog Number: 109676

## Product name

CYFIP2 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CYFIP2 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue, mouse brain tissue. Positive IHC detected in human gliomas tissue, human brain tissue. Observed molecular weight by Western-blot: 130 kDa

## Preparation

This antibody was obtained by immunization of CYFIP2 recombinant protein (Accession Number: XM\_011534516). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109676(CYFIP2 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:109676(CYFIP2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:109676(CYFIP2

# CYFIP2 antibody



Catalog Number: 109676

---

antibody) at dilution of 1:50 (under 40x lens)

IP Result of anti-CYFIP2 (IP:Catalog No:109676, 4ug; Detection:Catalog No:109676 1:700) with mouse brain tissue lysate 4000ug.

