

SLC3A1 antibody



Catalog Number: 115376

Product name

SLC3A1 antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SLC3A1 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human small intestine tissue. Positive WB detected in HEK-293 cells, HeLa cells, mouse kidney tissue. Positive IP detected in mouse kidney tissue. Observed molecular weight by Western-blot: 85 kDa

Preparation

This antibody was obtained by immunization of SLC3A1 recombinant protein (Accession Number: NM_000341). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115376(SLC3A1 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human kidney using Catalog No:115376(SLC3A1 antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded human kidney using Catalog No:115376(SLC3A1

SLC3A1 antibody



Catalog Number: 115376

antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-SLC3A1 (IP:Catalog No:115376, 4ug; Detection:Catalog No:115376 1:1000) with mouse kidney tissue lysate 4000ug.