

# LGALS3BP antibody



Catalog Number: 107376

## Product name

LGALS3BP antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

## Specificity

Human; other species not tested.

## Antibody description

LGALS3BP Mouse Monoclonal antibody. Positive IHC detected in human breast cancer tissue, human skin cancer tissue. Positive WB detected in human blood tissue. Positive IP detected in human plasma tissue. Observed molecular weight by Western-blot: 90 kDa

## Preparation

This antibody was obtained by immunization of LGALS3BP recombinant protein (Accession Number: NM\_005567). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


## Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:107376(LGALS3BP Antibody) at dilution of 1:2000



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:107376(LGALS3BP Antibody) at dilution of 1:50 (under 40x lens)




# LGALS3BP antibody



Catalog Number: 107376

Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:107376(LGALS3BP Antibody) at dilution of 1:50 (under 10x lens)



IP Result of anti-LGALS3BP (IP:Catalog No:107376, 5ug; Detection:Catalog No:107376 1:1000) with human plasma tissue lysate 4000ug.