

Synaptotagmin-4 antibody



Catalog Number: 115776

Product name

Synaptotagmin-4 antibody

Immunogen

[Human SYT4 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Synaptotagmin-4 Rabbit Polyclonal antibody.
Positive WB detected in human brain tissue.
Positive IHC detected in human gliomas tissue.
Observed molecular weight by Western-blot: 48kd

Preparation

This antibody was obtained by immunization of Synaptotagmin-4 recombinant protein (Accession Number: NM_020783). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115776(SYT4 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:115776(SYT4 antibody) at dilution of 1:100 (under 10x lens)