

Anti-LAT Rabbit antibody



Catalog Number: 167086

Product name

Anti-LAT Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

LAT Rabbit polyclonal antibody

Preparation

Antigen: Synthesized peptide derived from LAT .
at AA range: 190-270

Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA
and 0.02% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, ELISA

Dilutions

WB: 1/500 - 1/2000

ELISA: 1/10000

Validations



Western blot analysis of 293 using LAT antibody.
Secondary antibody was diluted at 1:20000