

# GATC antibody



Catalog Number: 110886

## Product name

GATC antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

GATC Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 15 kDa

## Preparation

This antibody was obtained by immunization of GATC recombinant protein (Accession Number: NM\_176818). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, IP


## Dilutions

Recommended Dilution:


IP: 1:200-1:2000

IHC: 1:20-1:200


## Validations



Immunohistochemistry of paraffin-embedded human liver slide using Catalog No:110886(GATC Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human liver slide using Catalog No:110886(GATC



Antibody) at dilution of 1:50

IP Result of anti-GATC (IP:Catalog No:110886, 3ug; Detection:Catalog No:110886 1:500) with mouse liver tissue lysate 12500ug.