

# YIPF5 antibody



Catalog Number: 116886

## Product name

YIPF5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

YIPF5 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, human heart tissue. Observed molecular weight by Western blot: 28kd

## Preparation

This antibody was obtained by immunization of YIPF5 recombinant protein (Accession Number: NM\_001024947). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116886(YIPF5 antibody) at dilution of 1:500