

# Anti-Poliovirus Receptor/PVR Rabbit antibody



Catalog Number: 167886

## Product name

Anti-Poliovirus Receptor/PVR Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB, IHC-P

## Antibody description

Poliovirus Receptor/PVR Rabbit monoclonal antibody

## Dilutions

WB: 1/1000

IHC: 1/50

## Preparation

Antigen: Recombinant protein of human Poliovirus Receptor/PVR

## Validations

## Formulation


50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal



Western blot detection of Poliovirus Receptor/PVR in HeLa,A549,U2OS cell lysates using Poliovirus Receptor/PVR Rabbit mAb(1:1000 diluted).Predicted band size:45kDa.Observed band size:60-80kDa.